



A new shape to your body

The pluses you do not expect:

Take advantage of 4 technologies:

- Lipolaser and Kotz waves: different types of technology with delivery via the same plate actuator. Elits' exclusive aesthetic technique for non-invasive treatments for Slimming, Modeling and Body Toning. One creates a lipolytic effect by inducing heat and the other simulates physical activity useful for the drainage of metabolic waste produced by overheating of adipose tissues.
- **Resistive radiofrequency**: technology useful for the stimulation of the production of new collagen. The result is a long-lasting lifting and firming effect.
- **Pressotherapy**: aesthetic treatment based on the external pressure of the areas relating to lower and upper limbs and abdomen for a draining result. The action takes place through a suit made up of leggings and an abdominal band divided up into fourteen combined sectors, automatically operated according to pre-established procedures.

For an effective action and visible results from the first treatments. The device works alone and in complete safety, **leaving the operator free to dedicate himself to other treatments or clients.**





Power: 230 Vac/50-60 Hz

Fuses: 2x5 Amp T 250 V

User Interface: Touch Screen Display

Dimensions and Weight:

Command Module: 37x42x40 cm - 21 kg

Base: 37x42x85 cm - 25 kg (optional)

Display: LCD Color 8" 800x600 pixels

Type of Apparatus:

BF Class

USB: Only for software/firmware

updates

Lipolaser	
Lipolaser wavelength	635 / 760 nm
Single plate maximum power	100 mW
Single led	18 mW/cm² max
Total plate - 10 LEDs	180 mW/cm² max
Applicators	6 plates
Kotz Waves	
Carrier frequency	2500 Hz
Maximum output voltage	120 Vac p-p +/- 10%
Maximum output current	1,1 mA/cm²
Carrier modulation	50 Hz
Duty cycle modulation	50% (10 ms ON 10 ms OFF)
T-ON	15 s
T-OFF	30 s (with burst 1 sec ON 2 sec OFF)
Applicators	6 plates
Resistive Multipolar Radiofrequency	
Bipolar plate output frequency	455 kHz +/- 10%
Maximum output voltage	180 Vac p-p
Frequency	455 kHz +/- 10%
T-ON	2 s
T-OFF	0,5 s
RF power	25 watt
Applicators	6 plates
Pressotherapy	
Max Pressure	60 mmHg per zone (if required, it can be increased up to 120 mmHg max)
	60 mmHg per zone (if required, it can be increased up to 120 mmHg max) 40 mmHg lower limbs - 20 mmHg abdomen - 30 mmHg upper limbs (if provided)
Max Pressure	
Max Pressure Standard Pressure	40 mmHg lower limbs - 20 mmHg abdomen - 30 mmHg upper limbs (if provided) Pressure bands for lower limbs and an abdominal band (if required, pressure
Max Pressure Standard Pressure Applicators	40 mmHg lower limbs - 20 mmHg abdomen - 30 mmHg upper limbs (if provided) Pressure bands for lower limbs and an abdominal band (if required, pressure bands for upper limbs)













Via E. Che Guevara, 26/C - 60022 CASTELFIDARDO (AN) ITALY

Email: info@elitsgroup.it / Tel: +39 071 720 67 39 elitsgroup.it